


Data file: D:\CHEM32\1\DATA\M7526Final.D
Sample name: M7526-Linacotide-02
Description: Mobile phase: Buffer A: 0.1% TFA in H₂O
Buffer B: 0.1% TFA in ACN
Gradient: 19-29% B in 20 min, Flow rate: 1 ml/min
Column: Unitary C18, 4.6*150 mm, 5 µm 100 Å
Instrument: Agilent HPLC 1100

Instrument: 1100

Location: 32
Injection: 1 of 1
Injection volume: 10.000
Acq. operator: SYSTEM



Signal: DAD1 A, Sig=220,4 Ref=360,100

RT[min]	Type	Width[min]	Area	Height	Area%	Name
8.391	VV B	0.2981	31.3064	1.4607	0.2597	
9.187	VV B	0.1870	11971.0615	976.5758	99.3100	
9.898	VV B	0.1071	13.9593	1.8708	0.1158	
10.084	VV B	0.1226	11.2764	1.2339	0.0935	
11.741	BV R	0.2189	26.6354	1.7371	0.2210	
Sum			12054.2389			