



Agilent Technologies

Sample name: M7546-SMER 28-01

Instrument: A010

Inj. volume: 1.5

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1D,Sig=230,4 Ref=off

RT [min]	Type	Width [min]	Area	Height	Area%
5.323	MM m	0.10	1.33	0.50	0.02
5.422	MM m	0.09	2.22	0.84	0.03
5.558	BM m	0.30	8512.94	2563.13	99.51
5.908	MM m	0.10	0.50	0.17	0.01
6.066	BM m	0.18	34.74	9.72	0.41
8.833	MM m	0.12	3.24	0.90	0.04