



Agilent Technologies

Sample name: M7727-CH-223191-04

Instrument: A010

Inj. volume: 2

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1B,Sig=210.0,4.0 Ref=off

RT [min]	Area	Height	Area%
7.422	19.54	6.10	0.17
8.692	4.16	1.33	0.04
10.470	7.38	1.64	0.06
10.711	3.87	1.02	0.03
10.865	5.74	1.48	0.05
11.791	11313.91	1535.28	99.64