



Agilent Technologies

Sample name: M7925-MST-312-02

Instrument: A012

Inj. volume: 2

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1B,Sig=210.0,4.0 Ref=off

RT [min]	Area	Height	Area%
2.069	16.33	5.31	0.12
2.318	35.70	11.68	0.27
2.533	5.69	1.37	0.04
7.883	13183.69	2520.28	99.20
8.417	44.06	10.93	0.33
9.193	5.20	1.44	0.04