



Agilent Technologies

Sample name: M7952-2-NP-01

Instrument: 1

Inj. volume: HPLC

Acq. method:

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1B,Sig=210,4 Ref=off

RT [min]	Area	Height	Area%
4.204	18.06	6.51	0.17
6.352	10578.11	2304.23	99.68
7.156	4.48	1.41	0.04
7.953	4.30	1.31	0.04
8.299	7.47	2.33	0.07