




Sample name: M8129-RN-1747-01  
Instrument: A010  
Inj. volume: 3  
Acq. method: HPLC  
Processing method:  
Manually modified:

Injection date:  
Location:  
Type:  
Sample amount:



Signal: DAD1B,Sig=210.0,4.0 Ref=off

RT [min]	Area	Height	Area%
5.640	3.96	1.34	0.03
6.120	5.36	1.07	0.04
6.269	37.19	9.79	0.29
6.339	8.40	3.20	0.06
6.460	6.60	2.22	0.05
6.533	6.21	1.97	0.05
6.826	2.19	0.71	0.02
7.018	33.90	11.70	0.26
7.100	3.32	1.22	0.03
7.367	12896.63	2600.27	99.18