


File D:\DATA\

Sample Name : M8254-SB297006-01

Method Info : Mobile Info:Mobile Phase:A:water(0.01%TFA) B:ACN(0.01%TFA)
Gradient: 5%B for 0.2min,increase to 95%B within 1.5min,95%B for 1.5min,back
to 5% within 0.01min
Flow Rate: 2ml/min
Column: SunFire ,50*4.6mm,3.5um
Temperature: 50 C



Integration Results for DAD1 A,Sig=214.4,Ref=off

RetTim	Width	Area	Height	Area%
1.98	0.01	696.13	789.54	100.00

Integration Results for DAD1 B, Sig=254.4 Ref=off

RetTim	Width	Area	Height	Area%
1.98	0.01	393.44	459.69	100.00

Ret. Time: 1.98

<<<< POSITIVE SPECTRA >>>>

