



Agilent Technologies

Project Name:

Operator:

Injection date:

Location:

Type:

Sample amount:

Sample name: M8891-AI-1-02


Instrument: A012

Inj. volume: 3

Acq. method: HPLC

Processing method:

Manually modified:



Signal: DAD1D,Sig=230.0,4.0 Ref=360.0,100.0

RT [min]	Type	Width [min]	Area	Height	Area%
6.800	MB m	0.25	198.32	57.19	1.06
7.137	MB m	0.14	18.61	5.72	0.10
8.051	MB m	0.13	12.81	3.92	0.07
8.158	BB	0.12	8.36	2.71	0.04
8.529	MB m	0.13	4.08	1.15	0.02
8.889	BM m	0.53	18494.96	2689.24	98.53
9.120	MM m	0.12	15.01	3.33	0.08
11.954	MM m	0.17	19.63	4.31	0.10