




Sample name: M8928-CEP-40783-01
Instrument: A012
Inj. volume: 1.5
Acq. method: HPLC
Processing method:
Manually modified:

Injection date:
Location:
Type:
Sample amount:



Signal: DAD1A,Sig=220,4 Ref=off

RT [min]	Peak Resolution USP	Area	Height	Area%
5.598		12.5	4.1	0.151
6.363		1.0	0.3	0.012
6.677		4.3	1.1	0.051
6.976		0.8	0.3	0.009
7.107		6.5	1.4	0.078
7.269		2.2	0.6	0.027
7.429		1.0	0.3	0.012
7.811		4.6	1.2	0.055
8.079		8268.5	1873.0	99.563
8.308		3.5	0.9	0.042