



Agilent Technologies

Sample name: M9016-Ezutromid-02

Instrument: A010

Inj. volume: 3

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1D,Sig=230.0,4.0 Ref=360.0,100.0

RT [min]	Type	Width [min]	Area	Height	Area%
4.942	BB	0.31	105.73	30.09	0.50
6.919	BB	0.14	4.62	1.21	0.02
8.159	BB	0.20	173.68	51.56	0.82
8.624	MB m	0.19	63.51	17.44	0.30
10.537	MB m	0.22	14.78	3.03	0.07
11.337	BV	0.83	20706.87	2381.79	98.13
11.658	VB	0.14	10.12	2.05	0.05
13.240	MB m	0.19	12.76	2.32	0.06
13.397	BB	0.15	6.11	1.45	0.03
13.675	MM m	0.13	3.78	0.85	0.02