



Agilent Technologies

Sample name: M9051-FCCP-03

Instrument: A010

Inj. volume: 2

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1A, Sig=220.0, 4.0 Ref=off

RT [min]	Peak Resolution USP	Area	Height	Area%
8.193		3.5	1.1	0.080
8.687		1.2	0.4	0.028
8.824		1.4	0.5	0.032
8.927		7.5	1.6	0.173
10.577		4335.9	671.1	99.431
11.832		9.6	2.3	0.220
11.969		1.5	0.4	0.035