



Agilent Technologies

Sample name: M9096-H2DCFDA-08

Instrument: A010

Inj. volume: 2

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1D, Sig=230.0,4.0 Ref=360.0,100.0

RT [min]	Type	Width [min]	Area	Height	Area%
9.504	BB	0.18	23.26	5.89	0.17
10.302	MB m	0.15	7.84	2.04	0.06
11.085	BM m	0.78	13520.68	2051.06	99.72
11.482	MM m	0.13	1.72	0.45	0.01
12.070	MB m	0.15	4.89	1.11	0.04