




Sample Name : M9148-SRT2104-01
Acq. Operator : 11
Spec. Reported : MS Integration
Acq. Method : D:\METHODS\2-POS-MON-1.M
Analysis Method : D:\METHODS\2-POS-MON-1.M
Sample Info : Easy-Access Method: '2-POS-MON-1'
Method Info : Mobile Phase: A: water(0.01%TFA) B: ACN(0.01%TFA)
Gradient: 5%-95% B in 1.5min, 95%B for 1.5 min
Flow Rate: 2ml/min
Column: Sunfire, 4.6*50mm, 3.5um
Oven Temperature: 50 C
MASS Range: 120-1000



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
1.67	0.01	6.17	6.98	0.31
1.83	0.02	1985.69	1865.11	99.23
2.31	0.02	9.30	7.31	0.46

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.83	0.02	1358.79	1290.81	99.38
2.31	0.02	8.44	7.39	0.62

Ret. Time: 1.85 <<< POSITIVE SPECTRA >>>

