



Agilent Technologies

Sample name: M9209-Belumosudil (SLx-2119)-01

Injection date:

Instrument: A010

Location:

Inj. volume: 2.5


Type:

Acq. method: HPLC

Sample amount:

Processing method:

Manually modified:



Signal: DAD1A, Sig=220.0, 4.0 Ref=off

RT [min]	Peak	Area	Height	Area%
5.893		0.9	0.3	0.007
6.007		27.0	3.9	0.188
6.287		23.8	5.6	0.166
6.565		14291.5	2496.0	99.594
7.606		2.1	0.7	0.015
9.418		4.4	1.1	0.030