



Agilent Technologies

Sample name: M9216-VcMMAE-02

Instrument: A010

Inj. volume: 2

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1B, Sig=210.0, 4.0 Ref=off

RT [min]	Area	Height	Area%
8.169	29.08	4.78	0.28
8.534	10177.56	2084.64	98.45
9.157	12.46	1.12	0.12
9.483	109.60	11.31	1.06
9.977	9.19	1.05	0.09