



Sample name: M9267-Coelenterazine-03

Injection date:

Instrument: A010

Location:

Inj. volume: 4


Type:

Acq. method: HPLC

Sample amount:

Processing method:

Manually modified:



Signal: DAD1E,Sig=260,4 Ref=360,100

RT [min]	Type	Width [min]	Area	Height	Area%
4.980	BM m	0.18	4.62	1.17	0.05
5.629	BB	0.15	10.41	3.15	0.12
6.589	MB m	0.11	4.40	1.31	0.05
6.770	MM m	0.09	48.62	17.00	0.58
6.877	MB m	0.36	8266.10	1628.51	98.14
7.267	BM m	0.12	10.73	2.48	0.13
7.335	MM m	0.13	10.65	1.87	0.13
7.521	MM m	0.04	0.20	0.13	0.00
7.611	MB m	0.16	21.08	5.30	0.25
7.741	BB	0.13	6.31	1.60	0.07
7.894	MM m	0.11	1.85	0.54	0.02
9.510	MM m	0.19	37.97	8.81	0.45