


Sample Name : M9356-AMG-510 (Sotorasib)-10 **Location** : P2-B-09
Acq. Operator : r **Inj** : 0
Spec. Reported : Integration **Inj Volume** : 2 ul



Integration Results for DAD1 A, Sig=214,4 Ref=off

| RetTim | Width | Area | Height | Area% |
|--------|-------|---------|---------|-------|
| 4.51 | 0.02 | 4.41 | 3.60 | 0.10 |
| 4.70 | 0.12 | 4302.99 | 2887.27 | 99.90 |