


=====

Acq. Operator : SYSTEM
Acq. Instrument : 6120
Sample Name : M9373-T-5224-03

Seq. Line : 3
Location : P2-A-05
Inj : 1
Inj Volume : 1.000 µl



=====

Area Percent Report

=====

Sorted By : Signal
Multiplier : 1.0000
Dilution : 1.0000
Do not use Multiplier & Dilution Factor with ISTDs

Signal 2: DAD1 B, Sig=254,4 Ref=off

Peak #	RetTime [min]	Type	Width [min]	Area [mAU*s]	Height [mAU]	Area %
1	7.823	MM	0.1010	6.76368	1.11660	0.2702
2	8.734	BB	0.0589	2496.23853	648.77179	99.7298

Totals : 2503.00220 649.88839