


MASS SPECTROMETRY REPORT



A: 3270.96±0.32

Sample Description

y
Analyst: SHAO
Sample: M9389-740 Y-P (PDGFR 740Y-P)-06
M.W.: 3270.7
Lot. No.:

Instrument

Probe:
Nebulizer Gas Flow:
CDL:
CDL Temp.:
Block Temp.:

Waters 2695

ESI
1.5L/min
-20.0v
250 °C
200 °C

Probe Bias: +4.5kv
Detector: 1.5kv
T. Flow: 0.5ml/min
B. Conc.: 50%H2O/50%ACN