



Sample name: M9490-YL-109-01

Instrument: A010

Inj. volume: 3

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1A, Sig=220.0, 4.0 Ref=off

RT [min]	Peak	Area	Height	Area%
7.081		28.8	5.5	0.145
7.929		27.9	5.7	0.140
8.241		24.1	5.4	0.121
8.347		0.0	0.0	0.000
8.906		286.6	74.7	1.443
9.259		19488.0	2851.6	98.150