

File ..DATA\T28\MONITOR\SVN-257-L-75144-LCMSA030.D Tgt Mass (EZX): 557.00
Seq. Line : 0
Sample Name : M9511-D0-264-01 Location : P1-A-09
Acq. Operator : T28-Monitor Inj : 1
Spec. Reported : MS Integration Inj Volume : 1.5 ul
Acq. Method : D:\METHODS\2-POS-15MIN.M
Analysis Method : D:\METHODS\2-POS-15MIN.M
Sample Info : Easy-Access Method: '2-POS-15MIN' 557.00
Method Info : Mobile Phase: A: water(10mM Ammonium hydrogen carbonate) B: ACN
Gradient: 0min 5%B,3min 5%B,10min 95%B,15min 95% B
Flow Rate: 1ml/min
Column: XBridge C18,4.6*150mm,3.5um
Oven Temperature: 45 C
MASS Range:110-1000



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
10.66	0.04	2098.93	748.87	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
2.97	0.01	4.50	9.08	0.29
3.89	0.01	5.16	8.67	0.33
3.99	0.01	7.22	9.33	0.46
10.66	0.04	1558.10	568.15	98.93

Ret. Time: 10.66

<<< POSITIVE SPECTRA >>>

