



Sample name: M9545-Hydrocortisone Hemisuccinate-01

Injection date:

Instrument: A012

Location:

Inj. volume: 3


Type:

Acq. method: HPLC

Sample amount:

Processing method:

Manually modified:



Signal: DAD1H,Sig=240,4 Ref=off

RT [min]	Type	Width [min]	Area	Height	Area%
6.520	MM m	0.15	146.99	41.84	1.28
7.197	MB m	0.33	11271.83	2826.45	98.39
7.635	MB m	0.12	4.96	1.23	0.04
7.789	MB m	0.14	9.41	2.12	0.08
7.967	MB m	0.09	2.78	0.97	0.02
8.058	BM m	0.09	2.50	0.89	0.02
8.629	MM m	0.14	3.83	0.90	0.03
9.049	MB m	0.14	5.69	1.46	0.05
9.316	MM m	0.16	8.21	1.67	0.07