


MASS SPECTROMETRY REPORT

M9555-Endothelin-1 (human, porcine)-04



A: 2491.29±0.29

Sample Description

Analyst: SHAO
Sample: M9555-Endothelin-1 (human, porcine)-04
M.W.: 2491.9
Lot. No.:

Instrument

Probe:
Nebulizer Gas Flow:
CDL:
CDL Temp.:
Block Temp.:

Waters 2695

ESI
1.5L/min
-20.0v
250 °C
200 °C

Probe Bias: +4.5kv
Detector: 1.5kv
T. Flow: 0.5ml/min
B. Conc.: 50%H2O/50%ACN