

Data file: D:\CHEM32\1\DATA\M9555-Endothelin-1 (human, porcine)-04.D
Sample name: M9555-Endothelin-1 (human, porcine)-04
Description: Mobile phase: Buffer A: 0.1% TFA in H₂O
Buffer B: 0.1% TFA in ACN
Gradient: 30-40% B in 20min, Flow rate: 1.0ml/min
Column: Unitary C18, 4.6*150mm, 5um 100A
Instrument: Agilent HPLC 1100


Instrument: 1100

Location: 22

Injection: 1 of 1

Injection volume: 10.000

Acq.operator: SYSTEM



Signal: DAD1 A, Sig=220, 4 Ref=360, 100

RT[min]	Type	Width[min]	Area	Height	Area%	Name
7.114	BB	0.2486	8.4403	0.5114	0.0645	
8.361	BV E	0.4838	148.7614	3.9824	1.1361	
9.031	VB R	0.3408	12936.8896	596.9920	98.7994	
Sum			13094.0913			