




Agilent Technologies

Sample name: M9605-Molnupiravir (EIDD-2801)-02 Injection date:
Instrument: A010 Location:
Inj. volume: 1 Type:
Acq. method: HPLC Sample amount:
Processing method:
Manually modified:



Signal: DAD1D,Sig=230.0,4.0 Ref=360.0,100.0

RT [min]	Type	Width [min]	Area	Height	Area%
3.497	BB	0.21	15.40	3.65	0.15
4.943	BB	0.65	10102.99	2745.89	99.85