



Agilent Technologies

Sample name: M9653-JMS-17-2-02

Instrument: A012

Inj. volume: 3

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1D,Sig=230.0,4.0 Ref=360.0,100.0

RT [min]	Type	Width [min]	Area	Height	Area%
6.986	BB	0.17	11.25	2.97	0.10
8.026	BV	0.44	10782.90	2600.34	99.85
8.339	VB	0.11	4.73	1.48	0.04