




Agilent Technologies

Sample name: M9727-SBI-115-01
Instrument: 1.5
Inj. volume: HPLC
Acq. method:
Processing method:
Manually modified:

Injection date:
Location:
Type:
Sample amount:



Signal: DAD1B, Sig=210,4 Ref=off

RT [min]	Area	Height	Area%
7.010	77.13	18.94	0.93
7.333	4.89	1.59	0.06
7.817	3.16	1.05	0.04
8.047	6.82	2.12	0.08
9.037	8177.00	1379.29	98.89
Sum	8268.99		