




Agilent Technologies

Sample name: M9739-GW-870086-01  
Instrument: 2  
Inj. volume: HPLC  
Acq. method:  
Processing method:  
Manually modified:

Injection date:  
Location:  
Type:  
Sample amount:



Signal: DAD1D,Sig=235,4 Ref=off

RT [min]	Type	Width [min]	Area	Height	Area%	Name
8.389	MM m	0.06	6.16	1.64	0.13	
12.851	MB m	0.11	4605.96	556.54	99.87	
	Sum		4612.12			