



Agilent Technologies

Sample name: M9938-Lazertinib-02

Instrument: A010

Inj. volume: 4

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1B, Sig=210.0, 4.0 Ref=off

RT [min]	Area	Height	Area%
7.481	6.46	2.58	0.08
7.612	7801.55	2089.68	99.64
7.837	4.56	1.52	0.06
8.381	8.03	2.13	0.10
9.055	9.34	2.46	0.12