




Sample name: M9957-Finerenone-01
Instrument: A010
Inj. volume: 2
Acq. method: HPLC
Processing method:
Manually modified:

Injection date:
Location:
Type:
Sample amount:



Signal: DAD1C,Sig=254.0,4.0 Ref=off

RT [min]	Peak	Width [min]	Area	Height	Area%
5.213		0.11	24.3	8.8	0.20
5.304		0.15	10.7	3.3	0.09
5.495		0.30	12099.4	3107.3	99.50
5.748		0.16	24.2	6.2	0.20
6.092		0.12	2.1	0.7	0.02