




9.493
9.454
8.415
8.402
8.388
7.366
7.360
7.328
7.321
7.314
7.307
7.300
7.260
7.254
7.245
7.234
7.224
7.220
7.191
7.188
7.180
7.176
5.040
4.352
4.345
4.340
4.337
4.332
3.898
3.884
3.872
3.857
3.171
3.164
3.148
3.141
2.802
2.778
2.766
2.749
2.742
2.725
1.915
1.904
1.843
1.832
0.825
0.816
0.805
0.708
0.696
0.684



Current Data Parameters
NAME MDL 28170-02
EXPNO 1
PROCNO 1

F2 - Acquisition Parameters
Date_ 20200613
Time 15.12 h
INSTRUM spect
PROBHD Z114607_0286 (
PULPROG zg30
TD 65536
SOLVENT DMSO
NS 16
DS 2
SWH 12019.230 Hz
FIDRES 0.366798 Hz
AQ 2.7263477 sec
RG 108.51
DW 41.600 usec
DE 6.50 usec
TE 298.0 K
D1 1.00000000 sec
TDO 1
SF01 600.1637060 MHz
NUC1 1H
P1 9.00 usec
PLM1 25.00000000 W

F2 - Processing parameters
SI 65536
SF 600.1600041 MHz
WDW EM
SSB 0
LB 0.30 Hz
GB 0
PC 1.00

0.94
0.94
11.00
2.03
0.97
1.01
0.98
1.02
1.00
3.18
3.02