

File D:\DATA\PRO2400-LCMS

Injection Date : 6 Nov 14 9:45 am +0800

Tgt Mass(EZX) :

Sample Name : OSI-930-LCMS

Location : P1-D-07

Acq. Operator : LCMS

Inj : 1

Spec. Reported : MS Integration

Inj Volume : 2 ul

Acq. Method : D:\METHODS\1-POS-3MIN.M

Analysis Method : D:\METHODS\1-POS-3MIN.M

Sample Info : *** No target masses specified ***

Method Info : Mobile Phase: A: water(0.01%TFA) B: ACN(0.01%TFA)

Gradient: 5%-95% B in 1.2min


Flow Rate: 2.2 ml/min

Column: Poroshell 120 EC-C18, 4.6*30mm, 2.7um


Oven Temperature: 50 C

Mass Range: 110-1500


*ADC1 A, ELSD




*DAD1 A, Sig=214,4 Ref=off



*DAD1 B, Sig=254,4 Ref=off



MSD1 TIC, MS File ES-API, Pos, Scan, Frag: 70



Integration Results for ADC1 A, ELSD

RetTim	Width	Area	Height	Area%
1.58	0.02	15.31	10.03	100.00

Integration Results for DAD1 A, Sig=214,4 Ref=off


RetTim	Width	Area	Height	Area%
1.59	0.02	944.26	568.08	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off


RetTim	Width	Area	Height	Area%
1.59	0.02	309.30	187.38	100.00

Ret. Time: 0.24


<<<< POSITIVE SPECTRA >>>>




Ret. Time: 0.63




Ret. Time: 1.08




Ret. Time: 1.12




Ret. Time: 1.59



Ret. Time: 2.18



Ret. Time: 2.26



Ret. Time: 2.32

