




Current Data Parameters
NAME Plerixafor-02
EXPNO 1
PROCNO 1

F2 - Acquisition Parameters
Date_ 20181030
Time_ 9.33
INSTRUM spect
PROBHD 5 mm PADUL 13C
PULPROG zg30
TD 65536
SOLVENT D2O
NS 8
DS 2
SWH 8012.820 Hz
FIDRES 0.122266 Hz
AQ 4.0894465 sec
RG 199.46
DW 62.400 usec
DE 6.50 usec
TE 296.3 K
D1 1.00000000 sec
TD0 1

==== CHANNEL f1 =====
SFO1 400.1324710 MHz
NUC1 1H
P1 10.00 usec
PLW1 20.00000000 W

F2 - Processing parameters
SI 65536
SF 400.1300000 MHz
WDW EM
SSB 0
LB 0.30 Hz
GB 0
PC 1.00



7.289

4.704

3.478
2.668
2.655
2.640
2.626
2.613
2.562
2.530
2.516
2.503
2.385
1.739
1.534

1.02

1.04

3.12

4.02

1.06

1.00

1.02

0 ppm