


File D:\DATA\  
Injection Date : 5 Jul 17 2:25 pm +0800 Tgt Mass (EZX) :  
Sample Name : R547-01 Location : P1-F-06  
Acq. Operator : LCMS Inj : 1  
Spec. Reported : MS Integration Inj Volume : 2 ul  
Acq. Method : D:\METHODS\1-POS-3MIN.M  
Analysis Method : D:\METHODS\1-POS-3MIN.M  
Sample Info : \*\*\* No target masses specified \*\*\*  
Method Info : Mobile Phase: A: water(0.01%TFA) B: ACN(0.01%TFA)  
Gradient: 5%-95% B in 1.2min  
Flow Rate: 2.2 ml/min  
Column: Poroshell 120 EC-C18, 4.6\*30mm, 2.7um  
Oven Temperature: 50 C  
Mass Range: 110-1500



## Integration Results for ADC1 A, ELSD

| RetTim | Width | Area | Height | Area%  |
|--------|-------|------|--------|--------|
| 1.38   | 0.02  | 2.59 | 2.09   | 100.00 |


## Integration Results for DAD1 A, Sig=214,4 Ref=off

| RetTim | Width | Area   | Height | Area%  |
|--------|-------|--------|--------|--------|
| 1.38   | 0.02  | 269.30 | 204.07 | 100.00 |


## Integration Results for DAD1 B, Sig=254,4 Ref=off

| RetTim | Width | Area   | Height | Area%  |
|--------|-------|--------|--------|--------|
| 1.38   | 0.02  | 308.93 | 225.76 | 100.00 |


Ret. Time: 1.34



Ret. Time: 1.39



Ret. Time: 1.83



Ret. Time: 1.91

