




LC/MS Report

Injection Date : 26 Nov 20
 Sample Name : RS 504393-02
 Method Info : Column:SunFire C18 (4.6x 50 mm,3.5 um)
 Mobile phase: H2O(0.01%TFA) (A) / ACN(0.01%TFA)(B)
 Elution program: Graient from 10 to95% of B in 1.4 min at 2ml/min
 Temperature: 50 °C
 Detection: UV (214 , 4 nm) and MS (ESI, Pos mode ,110 to 1000 amu)



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
1.52	0.02	1095.79	718.04	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.52	0.02	1011.45	673.80	100.00

Ret. Time: 1.52

<<<< POSITIVE SPECTRA >>>>

*ES-API Positive

