


Injection Date : 23 Feb 21
Sample Name : Rucaparib-01
Method Info Mobile Phase: A: water (0.01% TFA) B: ACN (0.01% TFA)
Gradient: 5% increase to 95%B within 1.3 min, 95%B for 1.7 min.
Flow Rate: 2.0ml/min
Column: Sunfire C18, 4.6*50mm, 3.5um
Column Temperature: 45 C
Detection: UV (214, 4nm) and MS (ESI, Pos mode, 110 to 1000 amu)



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
0.54	0.03	23.32	11.85	1.88
1.27	0.02	1218.87	936.60	98.12

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.27	0.02	674.04	521.67	100.00

Ret. Time: 1.27

<<<< POSITIVE SPECTRA >>>>

ES-API Positive

