

LCMS Report

Sample Information

LCMS


Acquired by : System Administrator
Sample Name : Sodium Phenylbutyrate-03
Method File : B 5%-95%-P-new - 2.6min.lcm

Date Acquired : 2021/6/29 10:47:40
Injection Volume : 5
Vial : 52


Column : Waters SunFire C18 50*4.6mm 5um 2.000ml/min 2.6min
Pump A : 0.1% FA and 10% ACN in H2O
Pump B : 0.1% FA and 10% H2O in ACN

Chromatogram


mAU



mAU



段号1



Integration Results for DAD

PDA Ch1 254nm

Peak#	Ret. Time	Height	Area	Area%
1	1.581	361869	606365	100.000
总计		361869	606365	100.000


PDA Ch2 214nm

Peak#	Ret. Time	Height	Area	Area%
1	1.582	2133828	5742904	100.000


Peak#	Ret. Time	Height	Area	Area%
总计		2133828	5742904	100.000

Mass Spectrum


Rt=0.287min 负




Rt=1.182min 正




Rt=1.429min 正




Rt=1.547min 负




Rt=1.542min 正




Rt=2.082min 正




Rt=2.132min 负



Rt=2.245min 负



Rt=2.402min 负



Rt=2.419min 正

