

File D:\DATA\PRO2772-LCMS

Injection Date : 12 Nov 14 10:52 am +0800

Tgt Mass (EZ) :

Sample Name : Sulfadimethoxine-LCMS

Location : P1-A-09

Acq. Operator : LCMS

Inj : 1

Spec. Reported : MS Integration

Inj Volume : 2 ul

Acq. Method : C:\Chem32\1\METHODS\1-POS-3MIN.M

Analysis Method : C:\Chem32\1\METHODS\1-POS-3MIN.M

Sample Info : Easy-Access Method: '1-POS-3min'

Method Info : Mobile Phase: A: water(0.01%TFA) B: ACN(0.01%TFA)

Gradient: 5%-95% B in 1.7 min


Flow Rate: 2.2 ml/min

Column: SunFire C18, 4.6*50mm, 3.5um


Column Temperature: 50 C

MASS Range: 110-900


*ADC1 A, ELSD




*DAD1 A, Sig=214,4 Ref=off



*DAD1 B, Sig=254,4 Ref=off



MSD1 TIC, MS File ES-API, Pos, Scan, Frag: 70



Integration Results for ADC1 A, ELSD

RetTim	Width	Area	Height	Area%
1.50	0.02	97.02	80.56	100.00

Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
1.50	0.02	4625.80	3124.79	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.50	0.01	2145.71	3065.18	51.03
1.51	0.01	2059.42	2969.18	48.97

Ret. Time: 1.52

<<<< POSITIVE SPECTRA >>>>

