




Injection Date : 22 Oct 20
Sample Name : Tanzisertib-04
Method Info : Column:SunFire C18 (4.6x 50 mm,3.5 um)
Mobile phase: H2O(0.01%TFA) (A) / ACN(0.01%TFA) (B)
Elution program: Gradient from 10 to95% of B in 1.4 min at 2ml/min
Temperature: 50 °C
Detection: UV (214 , 4 nm) and MS (ESI, Pos mode ,110 to 1000 amu)



Integration Results for DAD1 A, Sig=214,4 Ref=off


RetTim	Width	Area	Height	Area%
1.35	0.01	830.61	983.08	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off


RetTim	Width	Area	Height	Area%
1.35	0.01	2318.83	2597.73	100.00

Ret. Time: 1.35


<<<< POSITIVE SPECTRA >>>>




Ret. Time: 2.25



Ret. Time: 2.39




Ret. Time: 2.76



Ret. Time: 2.91

*ES-API Positive



Ret. Time: 2.96

*ES-API Positive

