

File D:\DATA\Temsirrolimus-HPLC

Acq. Method : D: HPLC

Analysis Method : D: HPLC

Sample Info : Easy-Access Method:

Method Info : Mobile Phase: A: water(0.01%TFA) B: ACN (0.01%TFA)

Gradient

Time 0 8 15


B % 5 95 95

Flow Rate: 1.0mL/min


Oven temperature: 40C

Column: SUNFIRE C18, 150*4.6mm, 3.5um

DAD1 A, Sig=214,4 Ref=off



DAD1 B, Sig=254,4 Ref=off



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
10.12	0.11	610.87	88.69	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
10.13	0.11	964.42	138.09	100.00