


7.608  
7.602  
7.585  
7.580  
7.435  
7.414  
7.400  
7.398  
7.346  
7.326  
7.307  
7.275  
7.271  
7.263  
7.254  
7.226  
7.114  
7.109  
7.090  
7.081  
7.059  
6.663  
6.307

3.775  
3.402  
2.266  
1.687  
-0.000

Sample Name: Tipifarnib-03  
 Chemical Formula: C<sub>27</sub>H<sub>22</sub>Cl<sub>2</sub>N<sub>4</sub>O  
 Molecular Weight: 489.40



CAS: 192185-72-1  
 CDCl<sub>3</sub>, 400 MHz

