




LC/MS Report

Injection Date : 6 Apr 21
 Sample Name : Zanubrutinib-01
 Method Info : Mobile Phase: A: water(0.01%TFA) B: ACN(0.01%TFA)
 Gradient: 5%-95% B in 1.5min
 Flow Rate: 2.0ml/min
 Column: SunFire C18 3.5um 4.6*50mm
 Oven Temperature:50 C
 Mass Range: 110-1000 UV: 214nm;254nm



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
1.82	0.02	1296.68	1304.77	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.82	0.02	756.64	763.81	100.00

Ret. Time: 1.82

<<<< POSITIVE SPECTRA >>>>

ES-API Positive

