

File D:\DATA\ABMOLE\HPLC01

Acq. Method : D:1-HPLC.M

Analysis Method : D:1-HPLC.M

Sample Info : Easy-Access Method: '1-HPLC.M'


Method Info : Column:Xbridge C18 (4.6x150 mm, 3.5um)

Mobile phase: H2O(10 mmol NH4HCO3) (A) / ACN(B)


Elution program: Gradient from 5 to 95% of B in 9min at 1ml/min Temperature: 50°C

Detection: UV (214 , 4 nm) and MS (ESI, Neg mode , 70 to 900 amu)

*DAD1 A, Sig=214,4 Ref=off



*DAD1 B, Sig=254,4 Ref=off



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
10.10	0.04	1178.31	506.07	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
10.10	0.04	896.93	386.36	100.00