

File D:\DATA

Injection Date : 18 Sep 15 2:45

Sample Name : LDE225 Diphosphate

Acq. Method : D:\METHODS\1-POS-3MIN.M

Analysis Method : D:\METHODS\1-POS-3MIN.M

Sample Info : *** No target masses specified ***

Method Info : Mobile Phase: A: water(0.01%TFA) B: ACN(0.01%TFA)

Gradient: 5%-95% B in 1.2min


Flow Rate: 2.2 ml/min

Column: Poroshell 120 EC-C18,4.6*30mm,2.7um


Oven Temperature: 50 C

Mass Range: 110-1500


*ADC1 A, ELSD




*DAD1 A, Sig=214,4 Ref=off



*DAD1 B, Sig=254,4 Ref=off



MSD1 TIC, MS File ES-API, Pos, Scan, Frag: 70



Integration Results for ADC1 A, ELSD

RetTim	Width	Area	Height	Area%
1.60	0.02	15.15	11.85	100.00

Integration Results for DAD1 A, Sig=214,4 Ref=off


RetTim	Width	Area	Height	Area%
1.60	0.02	1255.20	1180.07	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off


RetTim	Width	Area	Height	Area%
1.60	0.02	757.95	718.28	100.00

Ret. Time: 0.81


<<<< POSITIVE SPECTRA >>>>




Ret. Time: 1.14




Ret. Time: 1.30



Ret. Time: 1.37



Ret. Time: 1.61



Ret. Time: 2.39

