

Report

File D:\DATA

Injection Date : 28-May-15, 22:45

Sample Name : LT-674-01

Acq. Operator : HPLC

Acq. Method : C:\Chem32\1\METHODS\2-HPLC-15MIN.M

Analysis Method : C:\Chem32\1\METHODS\2-HPLC-15MIN.M

Sample Info :

Method Info : Mobile Phase: A: water (0.01%TFA) B: ACN (0.01%TFA)


Gradient: 0min 5%B, 3min 5%B, 10min 95%B, 15min 95%B

Flow Rate: 1.2 ml/min


Column: Eclipse XDB-C18, 4.6*150mm, 5um

Oven Temperature: 40 C

DAD1 A, Sig=214,4 Ref=off



DAD1 B, Sig=254,4 Ref=off



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
7.75	0.03	6.97	3.34	0.16
8.04	0.03	36.58	16.87	0.87
8.27	0.03	4179.84	1816.21	98.97

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
8.04	0.03	11.07	4.86	1.00
8.27	0.03	1100.08	487.79	99.00