



File D:\DATA\

Acq. Method : D: HPLC
Analysis Method : D: HPLC
Sample Info : Easy-Access Method:
Method Info : Mobile Phase: A: Water(0.01%TFA) B: ACN(0.01%TFA)
Gradient: 5%B increase to 95%B within 8min,
95%B for 7min, back to 5%B within 0.01min.
Flow Rate: 1 mL/min
Column: SunFire C18, 4.6*150mm, 3.5um
Column Temperature: 40 C



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
6.99	0.04	136.25	52.69	1.97
7.20	0.04	6792.45	2502.03	98.03

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
7.20	0.04	1778.70	668.00	100.00